



POWERHEART® AED G3 Pro
automated external defibrillator

ACCESSORIES



Introducing Cardiac Science Powerheart AED G3 Pro Accessories.

Designed specifically for medical professionals, Cardiac Science's latest line of accessories offers a variety of functions that enhance the capabilities of the Powerheart AED G3 Pro. The new series of accessories caters to the sophisticated demands of professional users by supporting great field and clinical usage of the Powerheart AED G3 Pro, and by simplifying the data management and customization processes.

The ECG Monitoring Kit facilitates the ability of medical professionals to monitor patients by providing continuous ECG monitoring for cardiac patients. Designed for long trips, transportation between hospitals, or constant cardiac monitoring, the ECG Monitoring Kit displays a non-diagnostic ECG heart rhythm on the Powerheart AED G3 Pro.

The Powerheart AED G3 Pro Rechargeable Battery and Charger combine convenience with reliability to afford professional users who demand more frequent use of the Powerheart AED G3 Pro. The Battery Charger quickly recharges the Rechargeable Battery to its full capacity in just three hours, and accepts all standard IEC main power cables.

To facilitate the professional user's ability to customize AED parameters and manage and transfer data, Cardiac Science introduces its infrared (IR) communication cable using the latest wireless technology.

Collectively or individually, Powerheart AED G3 Pro accessories enhance the professional user's capabilities. To learn more about Cardiac Science AEDs and accessories, visit our website www.cardiacscience.com.au, or email ds@cardiacscience.com.au, or call (03) 9429 2666.



CARDIAC SCIENCE

ORDER INFORMATION

PART NUMBER DESCRIPTION

5111-001

ECG Monitoring Kit

Provides Lead II viewpoint for accurate and continuous cardiac monitoring. Used in conjunction with the Powerheart AED G3 Pro, the ECG Monitoring Kit displays patient's non-diagnostic ECG heart rhythm on the AED screen. The electrode color convention is in accordance with the AHA standard and AAMI specifications.



9144-001

Powerheart AED G3 Pro Rechargeable Battery

Enables the Powerheart AED G3 Pro to deliver more shocks over its lifetime (60 shocks minimum, 100 is typical) and increases the duration of the ECG display (3 hours minimum, 6 is typical). The 1 lb. 3 oz. battery does not require any maintenance, and has a battery life of 2.5 years or 300 battery charge-discharge cycles, whichever come first.



9044-001

Powerheart AED G3 Pro Rechargeable Battery Charger and Cable

The new Powerheart AED G3 Pro Battery Charger quickly recharges the Powerheart AED G3 Pro Rechargeable Battery to its full capacity in just three hours (4.5 hours if rechargeable battery is completely depleted). The charger also accepts all standard IEC main power cables.



162-0108-001

IR Communications Cable

The wireless communication cable allows user to modify AED operational parameters, upload rescue data, and transfer information from three feet (0.9 meters) away. The IR Communications Cable is compatible with USB ports via an adapter.



CARDIAC SCIENCE

Nasdaq: DFIB

For more information contact Cardiac Science at:

Cardiac Science Australia

5 Prince Patrick Street, Richmond, VIC, 3121, Australia
Phone: (03) 9429 2666
Fax: (03) 9429 2611

(03) 9429 2666

WWW.CARDIACSCIENCE.COM.AU

EMAIL: DS@CARDIACSCIENCE.COM.AU

